

# AGP Report: Continuous Glucose Monitoring

Very High

13%

Goal: <5%

High

25%

Goal: <25%

Time Above Range

38%

Goal: >180

Time in Range

53%

Goal: >70%

Time in Tight Range

35%

Goal: >50%

Low

4.9%

Goal: <1%

Very Low

4.1%

Goal: <1%

Time Below Range

9.0%

Goal: <70

250

180

140

70

54

mg/dL

Each 1% time in range = ~15 minutes

Multiple daily injection (MDI) therapy

14 Days

Time CGM Active: 96.5%

Glucose Metrics

Average Glucose..... 160 mg/dL

Goal: <154 mg/dL

Glucose Management Indicator (GMI) ..... 7.1%

Goal: <7%


Glucose Variability ..... 44.6%

Defined as percent coefficient of variation

Goal: ≤36%

## Ambulatory Glucose Profile (AGP)

AGP is a summary of glucose values from the report period, with median (50%) and other percentiles shown as if they occurred in a single day.



## Daily Glucose Profiles

Each daily profile represents a midnight-to-midnight period.

